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Dear Sir, dear Madam, dear current or prospective LiMAX user

We are pleased to present today another edition of our monthly LiMAX Literature-Service. Each month we select one publication from international journals, which may be of interest to you regarding our unique LiMAX method.

This month we have selected the following publication:

Dosing of Antimycotic Treatment in Sepsis-Induced Liver Dysfunction by Functional Liver Testing with LiMAX

Kirchner et al, Case Reports in Critical Care, Dec. 2019

Critically ill patients with septic shock develop a liver dysfunction or liver failure in up to 50% of cases. Likewise, those patients often suffer from infections with Aspergillus or Candida, which require systemic antimycotic therapy. The predominantly hepatic metabolism of antimycotic agents evokes further hepatotoxic effects, which underlines the importance of correct dosing of these drugs. Hepatotoxicity is usually assessed by routine liver enzymes, offering, however, only static results. Earlier reports on intensive care therapy of septic patients have shown the capacity of the LiMAX test to detect liver dysfunction both early and effectively. In addition, a dosage control of antibiotics is possible. The presented case shows, for the first time, that an antimycotic therapy can be adapted patient specific based on the actual liver function.

The authors found the following results:

- A 56-year-old patient with pre-existing claudication was admitted with pain in the left leg due to an ischemia of the iliacal artery resulting from a newly detected infra renal aortic aneurysm
- After thrombectomy and subsequent vascular surgery the patient was transferred to the ICU in stable

condition

- On the 1st postoperative day, a bilateral renal infarction with an intra-aortic thrombus above the vascular surgical treatment and a compartment syndrome of the right leg developed, requiring both hemodialysis and fasciotomy
- On the 8th hospital day, a left sided colitis occurred, which was initially regressive but showed serious worsening on day 19, resulting in a septic shock with mechanical ventilation
- Following a re-laparotomy and Hartmann procedure the patient developed a multi organ failure with hemodynamic instability
- On day 25, liver enzymes were elevated, but the coagulation was normal; a LiMAx test revealed a good result of 351ug/kg/h, excluding a severe liver dysfunction
- Simultaneously a Candida sepsis evolved which required antimycotic Caspofungin therapy on day 28 with a loading dose of 70mg and a daily dose of 50mg
- A same-day LiMAx test indicated a liver dysfunction with a value of 201ug/kg/h
- 33 days after hospitalization another septic crisis and an increase of liver enzymes occurred, indicating a drug induced liver failure; consequently, the Caspofungin dose was reduced to 35mg
- An Aspergillus infection on day 41 was treated with continued Caspofungin therapy
- On the 48th day the liver function had improved as seen by a LiMAx value of 401ug/kg/h; the dose of Caspofungin was raised to 50mg again
- The antimycotic therapy was continued until day 65; on day 75 the patient was transferred to the normal ward without signs of infection and after 105 days in the hospital he was finally discharged

The authors conclude:

- Caspofungin is a potent antimycotic drug. However, to avoid overdosing and hepatotoxic effects the liver function should be monitored during such a therapy
- It is known that patients with Child-Pugh B classification should receive a reduced dose of Caspofungin, but how strong this reduction should be in such a patient with a score of 9 and existing hypalbuminemia is unknown
- The decision to further reduce the dose of Caspofungin to 35mg, justified by the decreased LiMAx value, was correct as seen by the quick recovery of the transaminases
- No other medication was reduced or interrupted during this period
- This case clearly shows the potential of a bedside monitoring of liver function in the systemic treatment of severe mycosis
- The LiMAx test detects a liver dysfunction several days prior to standard laboratory parameters and thus allows adjustment of drug dosages in real time

Humedics opinion:

- This is a very interesting case report from a clinical field in which the LiMAX test is still underrepresented
- On the other hand, critically ill patients in intensive care units present doctors with great challenges in assessing their organ functions
- For all major organs of a patient, there have long been real-time functional tests, only for the liver is this not the case.
- However, in the course of severe sepsis the liver is the organ whose failure or severe dysfunction can be decisive for survival, a timely and close monitoring of liver function is therefore necessary
- The LiMAX test is predestined for this purpose due to its uncomplicated availability and, at the same time, high explanatory power, practically presenting itself as the “ECG of the liver”

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: alexander.helmke@humedics.de), please unsubscribe [here](#).

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Humedics GmbH