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Dear Madam, Dear Sir, Dear current or prospective LiMAX user,

We are pleased to present the next edition of our quarterly LiMAX literature service today.

Each quarter, we select one or more publications from international journals, which may be of interest to you regarding our unique LiMAX method.

This time the following publication was selected:

Enzymatic liver function measured by LiMAX is superior to current standard methods in predicting transplant-free survival after TIPS implantation

Rashidi-Alavijeh et al., 2021. *Scientific Reports*, 11(1):13824, DOI: 10.1038/s41598-021-93392-5

Among patients with liver cirrhosis, portal hypertension (PH) is one of the main reasons for morbidity and mortality. In addition, with the occurrence of ascites and/or variceal bleeding (decompensated portal hypertension - DPH), mortality is rising.

One important treatment option for DPH is the placement of a transjugular intrahepatic portosystemic shunt (TIPS). This well-established therapeutic option for the reduction of portal hypertension complications is associated with different liver failure related complications such as hepatic encephalopathy (HE). Unfortunately, although several parameters such as age, MELD score, bilirubin levels, number of platelets and previous episodes of HE are known as risk factors for a TIPS placement, complications for this treatment still occur. Subsequently, there is still a need for prognostic markers for appropriate patient selection. Hence, the aim of this prospective study was to assess the predictive power of the LiMAX test on transplant-free survival of patients undergoing a TIPS placement.

In total, 30 patients suffering from decompensated cirrhotic portal hypertension and eligible for elective TIPS implantation were investigated at the Department of Gastroenterology and Hepatology at the University Hospital Essen between 2016 and 2018. In addition to angiographic measurements and TIPS placements, LiMAX tests were executed prior and after the TIPS implantation. Along with inferential statistical analysis, the areas under the receiver operation curves (AUROCs) were calculated.

The authors found the following results:

- Patients who survived without a liver transplantation exhibited higher median LiMAX values (1-day pre- and post-TIPS) and platelet counts as well as lower MELD scores than patients who died or received liver transplantation.
- LiMAX values decreased significantly after TIPS but the decrement resolved partly in 3-month survivors.
- While no correlation between LiMAX and hepatic venous pressure gradients (HVPG) values was found, a weak to moderate correlation between LiMAX values and MELD scores was identified. In addition, moderate correlations between LiMAX values and Child Pugh scores as well as the Freiburg index of post-TIPS survival (FIPS) occurred.
- Univariate analysis identified significant associations between LiMAX values (1-day pre and post TIPS) and 24-

week transplant-free survival.

- AUROC analysis showed good diagnostic accuracies for 1-day pre-TIPS LiMAx (sensitivity 100%, specificity 73.9%, AUROC 0.82, cut-off $\leq 205 \mu\text{g/kg/h}$), 1-day post-TIPS LiMAx (sensitivity 85.7%, specificity 78.3%, AUROC 0.85, cut-off $\leq 165 \mu\text{g/kg/h}$) and MELD values (sensitivity 71.4%, specificity 73.9%, AUROC 0.82, cut-off ≥ 15), while bilirubin as well as the post-TIPS LiMAx drop performed clearly worse.

The authors conclude:

- The impairment of liver function seems to be the main prediction factor for TIPS implantation related complications.
- The missing connection between LiMAx values and HVPG might be caused by a missing correlation between portal pressure and blood perfusion as well as subsequent neurohumoral and hemodynamic changes.
- As the LiMAx test showed significant associations with transplant-free survival, the LiMAx test could be beneficial for determining patient outcomes for TIPS candidates.
- While pre-interventional LiMAx values might constitute an adequate option for patient selection, post-interventional LiMAx values might help to prevent severe adverse events by exhibiting patients' needs for a more extensive and individual surveillance.
- Larger, multi-center studies are needed to confirm and verify study findings.

Opinion of Humedics:

- Once again it has been shown that the LiMAx test correlates with well-established as well as novel assessments such as the Child Pugh score or the novel Freiburg index of post-TIPS survival, respectively.
- We are particularly happy that in this study the LiMAx test has been shown to have a similar or even higher prediction on transplant-free survival than some of the routinely used parameters.
- Given that a reliable prognostic tool can support patient selection and reduce TIPS-associated adverse events, we believe that this study will help establish the LiMAx test as a prognostic tool in pre-TIPS placement.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAx Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: alexander.helmke@humedics.de), please unsubscribe [here](#).

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Humedics GmbH