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Dear Sir, dear Madam, dear current or prospective LiMAX user

We are pleased to present today our first edition of our new monthly LiMAX Literature-Service. Each month we select one or more publications from international journals, which may be of interest to you regarding our unique LiMAX method.

This month we selected the following publication:

Liver Maximum Capacity: A Novel Test to Accurately Diagnose Different Stages of Liver Fibrosis

Buechter et al, Digestion, October 2018

Chronic liver disease is – due to the increasing fraction of nutritional fatty liver – a major cause of mortality and morbidity worldwide and frequently lead to liver fibrosis followed by liver cirrhosis. As liver fibrosis is generally reversible, its early as possible diagnosis is of crucial importance for the further course of this disease. Liver biopsy as Gold Standard method due to its invasiveness and potential negative side effects is not suited as screening method, consequently alternative method based on biomarkers or elastography have been developed in the recent past.

Goal of this retrospective study was to investigate the accuracy of the LiMAX method with biomarkers and elastography to diagnose various stages of liver fibrosis.

In a total of 102 patients in whom a percutaneous or laparoscopic liver biopsy was taken, simultaneous LiMAX and elastography measurements were done and AAR, APRI, and FIB-4 scores determined.

The authors found the following:

- All non-invasive test showed a significant correlation to bioptic Desmet classification F1 – F4.

- For diagnosing the presence of fibrosis (Desmet ≥ 1) elastography showed highest sensitivity and specificity, immediately followed by LiMAX and FIB-4.
- In diagnosing significant (Desmet ≥ 2) or severe (Desmet ≥ 3) fibrosis, LiMAX was clearly superior to all other methods.
- Also, for the diagnosis of cirrhosis (Desmet = 4) LiMAX showed the highest sensitivity and specificity.

The authors conclude:

- LiMAX is superior to all established biomarker scores (AAR, APRI, FIB-4).
- For diagnosing the various stages of fibrosis LiMAX is suited better than elastography.
- To further increase sensitivity and specificity of a non-invasive diagnosis of cirrhosis compared to invasive biopsy a sequential algorithm of elastography and LiMAX is proposed.
- The combined use of both technologies may offer a significant and comfortable improvement in routine screening for liver disease.

Humedics opinion:

- After LiMAX in the field of liver surgery has already shown its supremacy to conventional methods, this is now one of the first publications to show that LiMAX will also in the field of Hepatology add significant benefit for monitoring and subsequent treatment of liver diseases.
- With the increasing incidence of digestive liver diseases (NASH) the importance of LiMAX, alone or in combination with elastography, will certainly grow further.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you want to unsubscribe from this Literature-Service (to: alexander.helmke@humedics.de) please click [here](#).

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Humedics GmbH