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Dear Sir, dear Madam, dear current or prospective LiMAx user

We are pleased to present today another edition of our monthly LiMAx Literature-Service.

Each month we select one publication from international journals, which may be of interest to you regarding our unique LiMAx method.

This month we selected the following publication:

Noninvasive Evaluation of Liver Function in Morbidly Obese Patients

Alizai et al, Gastroenterology Research and Practice, Feb 2019

The prevalence of overweight and obesity has dramatically increased over the last decades. Hepatic manifestation of morbid obesity is the non-alcoholic fatty liver disease (NAFLD) which can turn via non-alcoholic fatty liver inflammation (NASH) into Cirrhosis and Liver Cancer. Prevalence of NAFLD and NASH in morbidly obese patients is 50-90% and 10-50%, respectively, a screening for presence of those diseases therefore seems necessary. Common routine laboratory parameters are often normal in NASH, ultrasound has technical limitations in obese people, and other markers are so far not able to differentiate fatty liver from fatty liver inflammation. Biopsy is therefore still the method of choice for this differentiation, but because of invasiveness and potential complications not suitable as screening method. The authors investigated if the non-invasive LiMAX-Test is qualified for such differentiation. In 102 morbidly obese patients with mean Body-Mass-Index of 54kg/m² undergoing bariatric surgery, preoperatively a LiMAX test and NASH Clinical Scoring, intraoperatively a liver wedge resection followed by biopsy were performed and NAS (NAFLD Activity Score) computed.

The authors found the following results:

- 62% of patients showed LiMAX values lower than the normal value of 315ug/kg/h.
- 34% of patients showed no signs for NASH in biopsy, 39% showed borderline and 27% manifest signs for NASH (together 66%).
- LiMAX values significantly correlated negative with biopsically secured NAFLD Activity Score (NAS) and NASH Clinical Scoring according Campos et al.
- The LiMAX value in patients without NASH was significantly higher than in patients with borderline or manifest NASH.
- In patients with Diabetes type 2, however, LiMAX value was significantly lower than on non-diabetic patients.
- Correlations of LiMAX values to separate signs of liver cell damage (Ballooning, Inflammation, Steatosis, Fibrosis) were shown too.
- The Cut-Off value for differentiation of No-NASH vs. manifest or borderline NASH was 288ug/kg/h with sensitivity / specificity of 85% / 83% and AUROC of 0.859.

The authors conclude:

- Based on these results, LiMAX test seems to be qualified as screening method for presence of NASH in obese patients.
- LiMAX testing revealed high sensitivity and specificity to diagnose NASH in morbidly obese patients and results correlate well with histology results and clinical scoring systems.
- LiMAX thus can facilitate NASH screening in patients at risk, furthermore, repeated LiMAX tests allow continuous monitoring of the disease course or prove success of therapeutic interventions.

Humedics opinion:

- This publication shows clearly that LiMAX is not only suited for determining respectively differentiating the severity of, but also as screening method for detecting liver disease.
- For differentiation of a harmless stage of pure liver cell adiposis (NAFLD) and beginning or manifest inflammation (NASH) because of its robust Cut-Off values, rapid availability, and simple application LiMAX is superior to all other methods like NAFLD score, elastography, MRT, or biopsy.
- Biopsy and MRT of course allow an even better differentiation of the various

stages of liver cell damage, however they do tell nothing about current, actual liver function.

- Because especially the latter – not only in bariatric patients – is responsible for the further course of the disease and thus a patients' outcome, liver function should be regularly monitored by the LiMAX test.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you want to unsubscribe from this Literature-Service (to: alexander.helmke@humedics.de) please click [here](#).

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Humedics GmbH