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Dear Sir, dear Madam, dear current or prospective LiMAX user

We are pleased to present today another edition of our monthly LiMAX Literature-Service.

Each month we select one publication from international journals, which may be of interest to you regarding our unique LiMAX method.

This month we selected the following publication:

**The utility of the 13C-methacetin breath test in predicting the long-term survival of patients with decompensated cirrhosis**

*Moran, et al, J Breath Research 11 (2017)*

During the natural course of liver cirrhosis there are 2 phases with different survival rates; the compensated phase with a hepatic venous pressure gradient less than 10mmHg and the decompensated phase with a pressure gradient higher than 10mmHg and complications like variceal bleeding, ascites, and encephalopathy. The mainly used survival predictors for both phases are the Child-Pugh Score for long term and the MELD Score for short term survival. Both scores however are limited by individual judgement of the clinical parameters, inconsistency of the biochemical measurements and the fact that none of the scores includes functional liver reserve. Methacetin breath tests could already show in the past that they are correlated to the Child-Pugh score as well as with the severity of the liver disease. While the Child-Pugh Score is also very frequently used to estimate the long-term survival prognosis of patients with decompensated cirrhosis, there is little knowledge in this respect about the Methacetin breath test.

The authors therefore investigated the usefulness of a Methacetin breath test for survival prediction of patients with decompensated cirrhosis.

123 patients with decompensated liver cirrhosis were prospectively included in the study. Reasons for cirrhosis were alcohol abuse, chronic hepatitis B or C, auto immune hepatitis, primary biliary cholangitis, non-alcoholic fatty liver disease and onset cryptogenic liver disease. Upon inclusion in the study the MELD and Child-Pugh scores and routine liver laboratory parameters were determined and a Methacetin breath test was carried out; patients were followed-up for 3 years. The Methacetin breath test was performed using an experimental setup with an oral intake and the target parameter was the  $^{13}\text{C}$  increase after 15 minutes (DOB15).

***The authors found the following results:***

- Of the 123 enrolled patients 30 died for hepatic reasons, 7 died from non-hepatic reasons and were censored at the time of death, total follow-up time was 28.8 months.
- Mortality after 1 year was 5.8%, after 2 years 9.7% and 10.8% during the third year.
- The Cut-Off values for the 3 year mortality prediction were a Child-Pugh Score  $\geq 8$ , a MELD score  $\geq 12$  and DOB15  $\leq 4.5_{0/00}$ .
- Survival probability after 1 year was 0.89 for patients with DOB15  $\leq 4.5_{0/00}$  and 0.98 for patients with DOB15  $> 4.5_{0/00}$ , after 2 years 0.76 with DOB15  $\leq 4.5_{0/00}$  and 0.92 with DOB15  $> 4.5_{0/00}$ . After 3 years survival probability was 0.57 with DOB15  $\leq 4.5_{0/00}$  and 0.85 with DOB  $> 4.5_{0/00}$ .
- Survival curves were significantly different for the groups with DOB  $\leq 4.5_{0/00}$  and  $> 4.5_{0/00}$ .
- In a multivariate analysis mortality was 2.58 times higher in patients with DOB 15  $\leq 4.5_{0/00}$  compared to patients with DOB15  $> 4.5_{0/00}$ .

***The authors conclude:***

- The Methacetin breath test offers an accurate determination of the hepatic functional reserve and allows a survival prediction for patients with decompensated liver cirrhosis.
- The results of the Methacetin breath test were similar to those for Child-Pugh and MELD Scores, the Methacetin breath test is easy and quick to perform.
- The Methacetin breath test thus qualifies as a screening method in patients at risk of liver disease, as a follow-up method in liver cirrhosis and as a criterion for liver donors.
- Before the method can be integrated into routine medical practice, commercial and clinical challenges need to be mastered.

***Humedics opinion:***

- In light of the FDA changing its position that liver drug trials need to show histologic improvements as a clinical endpoint, this already 2 year old study receives a totally new significance.
- As this study shows, the Methacetin breath test has the potential to predict the long-term survival of patients with liver disease based on an initial test result. Should drug therapy of liver disease lead to an improvement in the result of a Methacetin breath test, there is a high likelihood that the long term survival of the patient improves as well.
- The Methacetin breath test technique used in this study was an experimental one which although similar to the LiMAX test, is technologically inferior.
- By using the LiMAX technology and especially intravenous, instead of error-prone oral, Methacetin administration, the results of this study would most probably be even clearer.
- The cut-off value of DOB15  $0.45_{0/00}$  corresponds approximately to a LiMAX value of 75 ug/kg/h, similarly low cut-off values for prediction of survival have already

been shown in the LiMAX literature many times.

- Since the LiMAX test has already been successfully used in hospitals for several years, the authors' request to overcome clinical and commercial hurdles can now be regarded as fulfilled.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: [alexander.helmke@humedics.de](mailto:alexander.helmke@humedics.de)), please unsubscribe [here](#).

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Humedics GmbH