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Dear Sir, dear Madam, dear current or prospective LiMAX user

We are pleased to present today another edition of our monthly LiMAX Literature-Service. Each month we select one publication from international journals, which may be of interest to you regarding our unique LiMAX method.

This month we have selected the following publication:

**Gd-EOB-DTPA-enhanced MRI T1 relaxometry as an imaging-based liver function test compared with 13C-methacetin breath test**

Theilig et al, Acta Radiologica, July 2017

Gd-EOB-DTPA (Gadolinium) enhanced magnetic resonance imaging (MRI) and especially MRI T1 Relaxometry is a procedure allowing the assessment of spatial distribution of liver function without ionizing radiation and can be easily integrated into preoperative clinical routine. This method has, in the past, been compared to various other methods for the assessment of liver function, including Child-Pugh classification, MELD Score, ICG-elimination rate, scintigraphy and direct Fibrosis staging. A comparison with the functional assessment of the actual liver function by the LiMAX test has not yet been performed and was therefore the aim of this study.

From 597 patients who underwent Gadolinium-enhanced MRI in 2015 and 2016, 53 patients were retrospectively identified in whom a LiMAX test was performed within a period of less than 30 days before or after the MRI. For these patients in addition to the LiMAX test and the bilirubin value, total liver volume, functional liver volume, and various T1 relaxometry parameters pre and post Gadolinium injection during the MRI were determined and compared to each other.

***The authors found the following results:***

- 31 patients had a LiMAX value higher than 315 ug/kg/h, 22 patients had a reduced value between 315 and 140 ug/kg/h and 8 patients had a LiMAX value lower than 140 ug/kg/h.
- All MRI T1 relaxometry parameters showed a significant correlation to the LiMAX value, bilirubin showed as expected no correlation.
- In ROC analyses for prediction of a LiMAX value higher or lower than 315 ug/kg/h the majority of T1 relaxometry parameters also showed satisfactory results, good results were also found for some of these parameters for the prediction of a LiMAX value lower than 140 ug/kg/h; T1 post and T1 RR were leading in all ROC analyses.
- In a multiple linear regression analysis, only the MRI T1 RR value was found to be a good predictor of LiMAX liver function, whereas total and functional liver volume as well as bilirubin and T1 post showed no significant results.

***The authors conclude:***

- All MRI T1 relaxometry parameters showed a correlation to the LiMAX value and allow a discrimination of LiMAX values higher than 315 ug/kg/h, of 315 ug/kg/h to 140 ug/kg/h and below 140 ug/kg/h.
- While the T1 RR value was best for the prediction of a good to sufficient LiMAX value higher than 140 ug/kg/h, the T1 post value was the best for the prediction of impaired liver function lower than 140 ug/kg/h.
- The unclear result for T1 RR and T1 post may be caused by the small number of cases of only 8 patients with a LiMAX value lower than 140 ug/kg/h.
- Neither the functional nor the total liver volume could contribute to the prediction of the LiMAX value.
- The comparison of the Gadolinium-enhanced MRI to the LiMAX test is superior to the more often performed comparison to the ICG elimination rate because Gadolinium and ICG are excreted in the liver by the same metabolic pathway, whereas the Methacetin of the LiMAX test is metabolized differently.
- It should be noted that in context of this study predominantly the excretory liver function has been investigated.
- If MRI Gadolinium T1 relaxometry proves to be equal to the LiMAX test, this may have an impact on preoperative procedures since MRI would also provide information about pathology and spatial distribution of liver function.
- Further studies are required to confirm these results, in particular the preoperative prediction of postoperative liver function using Gadolinium-enhanced MRI should be further investigated.

***Humedics opinion:***

- Basically, it is tempting to think that a preoperative MRI examination will not only gather information on the anatomy and pathology of the liver, but also on its function.
- However, before such a method could enter clinical routine many more studies are required, as already stated by the authors.
- Whether the results will be confirmed at all needs anyway first to be seen; the fact that Gadolinium and ICG are eliminated via the same metabolic pathway leads us to expect that in further comparisons of this kind the limitations of the

ICG elimination rate compared to the LiMAX test will also be found for MRI T1 relaxometry.

- Completely disregarded is the fact that – in contrast to an MRI, which is unpleasant for patients - the LiMAX test can be repeated, also postoperatively, in a simple manner and at any time.
- Also, Gadolinium is a not a completely harmless contrast agent as demonstrated by the recent warnings and limitations on its use issued by the European health authority. This MRI procedure is therefore not really suitable for pre- and postoperative monitoring.
- In individual cases before liver surgery it may therefore be conceivable to replace the LiMAX test with Gadolinium-enhanced MRI examinations, however in the general routine we consider more sensible for LiMAX and MRI to complement each other.
- Especially outside purely surgical questions, such as in Hepatology/Gastroenterology, but also regarding the question about whether surgery can and should be performed at all, the LiMAX test will continue to expand its standing because it is easy to perform and provides clear informative value and rather act as a gatekeeper before an MRI than be replaced by the same.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: [alexander.helmke@humedics.de](mailto:alexander.helmke@humedics.de)), please unsubscribe [here](#).

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Humedics GmbH