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Dear Sir, dear Madam, dear current or prospective LiMAX user

We are pleased to present today another edition of our monthly LiMAX Literature-Service. Each month we select one publication from international journals, which may be of interest to you regarding our unique LiMAX method.

This month we have selected the following publication:

Predictive Power of Liver Maximum Function Capacity Test in Transjugular Intrahepatic Portosystemic Shunt Patients: A Pilot Study

Reichert et al, Dig Diseases, Oct. 2019

A transjugular, intrahepatic shunt (TIPSS) is performed in patients with various advanced liver diseases, the procedure itself however often leads to complications associated with liver function, too. Hepatic encephalopathy which is related to increased morbidity and mortality, for example, occurs in approximately 60% of patients and is a drawback of this life prolonging intervention. Hepatic decompensation is less common but constitutes a life-threatening complication as well.

Although a variety of clinical parameters and scores (i.e. MELD and Child-Pugh) is available for stratification of patients at risk, liver function associated complications still occur after TIPSS. Various quantitative liver function tests have been investigated in the past as well, the LiMAX test however has not yet been evaluated in a TIPSS context.

The aim of this study was therefore to measure liver function with the LiMAX test before and after TIPSS and to prospectively assess potential risk factors for the development of liver function associated complications after TIPSS.

The authors found the following results:

- In 20 patients the LiMAX value was measured before and after TIPSS implantation, in 5 patients another LiMAX value was measured after 14 days
- Child-Pugh, MELD score and other laboratory parameters were determined before and after TIPSS as well, 90% of patients were graded as Child-Pugh class B
- 8 patients (40%) developed a liver function associated complication
- Patients who developed such a complication showed a significantly higher hepatovenous pressure gradient before the TIPSS procedure
- Patients suffering from liver function related complications had a significantly larger drop of the LiMAX value after TIPSS implantation compared to patients without complications

- In a multivariate analysis the drop of the LiMAX value was the only parameter which was significantly associated with such complications, neither bilirubin nor INR increase showed a significant relationship
- A drop of the LiMAX value by 58ug/kg/h as cut-off value can predict a liver function associated complication with a sensitivity of 0.75 and a specificity of 0.91, the bilirubin increase did not show any significant predictive power

The authors conclude:

- The LiMAX value did not correlate to the hepatovenous pressure gradient, consequently the LiMAX value is not influenced by portal hypertension
- A marked drop of the LiMAX value, however, predicts the imminent development of liver function related complications after a TIPSS procedure
- Neither the level, nor the increase of bilirubin predicted a potential complication
- The existing literature describes a relationship between a high hepatovenous pressure gradient or a strong decrease of the latter and liver function associated complications; this was the case in this study as well, however without statistic significance
- A high MELD score however was not associated with complications in this study. This may be due to the relatively low MELD scores of the enrolled patients, similar results were found for the Child-Pugh Score
- Higher age and previous hepatic encephalopathy have been described as risk factors as well, which could also not be confirmed in this study
- Prospective LiMAX measurements before and after TIPSS can detect a potential strong drop of the liver function after TIPSS
- LiMAX allows to predict the development of perinterventional complications and could be helpful to develop preventive strategies to avoid TIPSS complications

Humedics opinion:

- As **also shown recently for TACE procedures**, this study confirms that a routine use of the LiMAX test before and after interventional procedures can be helpful to prevent or reduce complications.
- Although case numbers in these studies are small, they clearly show the strong potential a LiMAX monitoring has for such interventions.
- We are convinced that future studies in this area with higher case numbers will confirm these first results and that the LiMAX test will be included in the guidelines for interventional liver therapies

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: alexander.helmke@humedics.de), please unsubscribe [here](#).

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Humedics GmbH