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Dear Sir, dear Madam, dear current or prospective LiMAX user

We are pleased to present today another edition of our monthly LiMAX Literature-Service. Each month we select one publication from international journals, which may be of interest to you regarding our unique LiMAX method.

This month we have selected the following publication:

### **Dynamic determination of functional liver capacity with the LiMAX test in post cardiac arrest patients undergoing targeted temperature management – a prospective trial**

***Nee et al., Acta Anaesthesiol Scand, April 2020***

An increase in transaminases, particularly of ASAT, is frequently seen in survivors of cardiac arrest. ASAT levels higher than 10.000 U/L define a hypoxic hepatitis as a secondary form of acute liver failure associated with shock and hypoxemia. The characteristics of hypoxic hepatitis are not well known yet, the incidence however is reported with app. 11 to 14 percent.

In comatose survivors of cardiac arrest, treatment often involves a targeted reduction in body temperature, but this has an impact on pharmacokinetics and pharmacodynamics resulting in elevated drug levels or reduced drug effects. The CYP450 1A2 enzyme system plays a central role in the hepatic drug metabolism and can be measured directly with the LiMAX test.

The aim of this pilot study was to investigate the functional liver capacity in relationship to the severity of hypoxic hepatitis in patients after cardiac arrest, treated with targeted lowering of their body temperature.

#### ***The authors found the following results:***

- A total of 32 patients with a mean age of 65 were prospectively enrolled in the study
- 5 patients received amiodarone (a weak CYP1A2 inhibitor) before, 4 during their ICU, 1 patient received ciprofloxacin (strong CYP 1A2 inhibitor)

- 18 patients were treated with insulin (a weak CYP 1A2 inducer), 3 patients had a history of tobacco use (potential CYP 1A2 induction)
- All other drugs used for sedation or analgesia are not known as inducer or inhibitor of CYP 1A2
- The liver enzymes ASAT, ALA, and GLDH were markedly elevated and normalized in the course of 10 days, bilirubin, coagulation factors and lactate were normal at all times, albumin was slightly lowered at all times
- LiMAx values on day 1 (at 33° C body temperature), on day 2 (during the warm-up phase) and on day 5 after successful resuscitation were below the normal range
- The minimum LiMAx value was seen on day 2, a recovery was seen on day 10, all changes of the values between measurements were statistically significant
- In the multiple regression analysis, a significant relationship of resuscitation length and the ASAT level on day 1 was seen
- No relationship was seen between the lowest LiMAx value on day 2 and the resuscitation time, nor was there a significant correlation between the resuscitation time and the differences in LiMAx values between day 1 and 2 and day 2 and 5 after resuscitation

***The authors conclude:***

- CYP 1A2 activity was significantly lowered at day 2 after resuscitation and slowly recovered until day 10
- The duration of resuscitation significantly impacted the static hepatic parameters, but not the functional parameters
- The results seen for liver enzymes ASAT and ALAT were similar to those known from the literature
- No patient in this cohort met the criteria for hypoxic hepatitis, but the pattern was typical of transient, mild hypoxic liver damage
- The duration of resuscitation-related hypoxia does not seem to be the main reason for the reduced liver function shown by the LiMAx test, it is probably more due to a combination of reperfusion-related systemic inflammation, severity of hypoxic liver damage and body temperature
- During the study period some patients received drugs that are known to interact with the activity of CYP 1A2, a potential impact seems to be unlikely, however cannot be ruled out completely
- The metabolism of sedatives and other drugs can be influenced by reduced liver function during the initial intensive care therapy; this should be taken into account in clinical decision-making when administering drugs
- Further clinical studies should investigate this effect and the influence of functional liver capacity on drug dosing to further optimize therapy following cardiac arrest

***Humedics opinion:***

- Earlier studies have already investigated and shown how the measurement of liver function with the LiMAx test can help optimize drug dosage and possibly improve patient outcome when treating critically ill patients in intensive care units
- This study is the first to show the impact of cardiac arrest and subsequent intensive care therapy using body temperature management on liver function
- Although it is too early to make concrete monitoring recommendations based on the LiMAx test in intensive care units, the interest of clinicians is obvious, and we very much look forward to further upcoming studies

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAx Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: [alexander.helmke@humedics.de](mailto:alexander.helmke@humedics.de)), please unsubscribe [here](#).

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Humedics GmbH