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Dear Madam, Dear Sir, Dear current or prospective LiMAx user,

We are pleased to present the next edition of our quarterly LiMAx literature service today.

Each quarter, we select one or more publications from international journals, which may be of interest to you regarding our unique LiMAx method.

This time the following publication was selected:

### **Volume–Function Analysis (LiMAx Test) in Patients with HCC and Cirrhosis Undergoing TACE – A Feasibility Study**

**Reichert et al., 2020. *Digestive Diseases and Sciences*, DOI: 10.1007/s10620-020-06535-5.**

Hepatocellular carcinoma (HCC) - the third most common cause for cancer-dependent death - ranges among the most frequent tumors worldwide. Its incidence differs between geographical regions and significantly increased within the last decade.

According to current international guidelines, transarterial chemoembolization (TACE) – also used for bridging patients to liver transplantation – is the gold standard for HCC in intermediate stage. As several parameters like Child-Pugh stages and serum bilirubin level have been identified as risk factors for post-TACE liver failures and as several ambitions have been undertaken to develop and evaluate scoring systems and algorithms for HCC patients, the aim of this prospective feasibility study was to assess the relationship of LiMAx values and other HCC-related parameters in HCC patients.

In total, 11 HCC patients with cirrhosis undergoing repetitive TACE at the Department of General, Visceral, Vascular and Pediatric Surgery Homburg were investigated. Etiology of cirrhosis, several laboratory measurements like bilirubin, albumin and sodium, Child-Pugh Scores (CPS), model for end-stage liver disease (MELD) scores and LiMAx tests were carried out 1 day pre and post TACE treatment. Univariate analyses were performed to investigate group differences.

#### **The authors found the following results:**

- Pre-TACE LiMAx values correlated with pre-TACE ALBI scores and albumin levels.
- While a non-significant trend occurred between pre-TACE LiMAx values and pre-TACE bilirubin, pre-TACE LiMAx values significantly correlated with post-TACE bilirubin levels.
- Serum bilirubin increased significantly after TACE. This increase was significantly stronger for pre-TACE LiMAx values  $\leq 150$   $\mu\text{g}/\text{kg}/\text{h}$  than for those  $> 150$   $\mu\text{g}/\text{kg}/\text{h}$ .
- A statistical trend towards slightly reduced LiMAx values after TACE occurred.
- Higher bilirubin as well as INR levels, larger spleen sizes and lower platelet counts were observed in patients with lower LiMAx ( $\leq 150$   $\mu\text{g}/\text{kg}/\text{h}$ ) values.

#### **The authors conclude:**

- Higher bilirubin as well as INR levels, larger spleen sizes and lower platelet counts confirm the well-known relation of reduced LiMAx levels and increasing stages of cirrhosis.
- LiMAx measurements are reliably reproducible in patients with HCC and cirrhosis undergoing multiple TACE

treatments.

- LiMAX levels did not correlate with ALBI levels, probably due to the similar stage of cirrhosis of the included patients.
- Even if the relationship of reduced LiMAX values and post-TACE bilirubin increments might reflect changes in vascularization after TACE, future studies are required to investigate if LiMAX values could be an important marker to predict post-TACE hepatic failure (PTHF).
- Despite some limitations like the small sample size and the absence of major hepatic failure, significant correlations with surrogate markers indicate that LiMAX is a robust method for liver function assessment in cirrhotic liver. Nevertheless, the optimal time point of LiMAX measurements after TACE has to be investigated in further studies.

### **Opinion of Humedics:**

- We are pleased that the study has shown first insights in the use of the LiMAX test in patients undergoing TACE. We are happy that the relationship between LiMAX test values and surrogate parameters indicates that the LiMAX test might have potential as predictor of PTHF.
- Similar to many other studies, this study exhibits that LiMAX values are strong predictors for several other liver function related laboratory parameters. It further underpins the approach to separate participants for different indications by predefined LiMAX thresholds.
- We would appreciate further research investigating the use of the LiMAX test for transarterial chemoembolization.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: [alexander.helmke@humedics.de](mailto:alexander.helmke@humedics.de)), please unsubscribe [here](#).

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Humedics GmbH